

Antibodies Product Datasheet

Product Name:	GPX8 Antibody
Catalog Number:	C48977
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPX8 AA range: 158-172
Concentration:	1 mg/ml
Gene Name:	GPX8 UNQ847/PRO1785
Protein Name:	Probable glutathione peroxidase 8 (GPx-8) (GSHPx-8) (EC 1.11.1.9)
Other Name:	
Human Swiss Prot No.:	Q8TED1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPX8 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year