

Antibodies Product Datasheet

Product Name:	GM2A Antibody
Catalog Number:	C48998
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GM2A
Concentration:	1 mg/ml
Gene Name:	GM2A
Protein Name:	Ganglioside GM2 activator (Cerebroside sulfate activator protein) (GM2-AP) (Shingolipid activator protein 3) (SAP-3) [Cleaved into: Ganglioside GM2 activator isoform short]
Other Name:	
Human Swiss Prot No.:	P17900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GM2A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year