

Antibodies Product Datasheet

Product Name:	IRF2BP1 Antibody
Catalog Number:	C49001
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IRF2BP1
Concentration:	1 mg/ml
Gene Name:	IRF2BP1
Protein Name:	Interferon regulatory factor 2-binding protein 1 (IRF-2-binding protein 1) (IRF-2BP1) (Probable E3 ubiquitin-protein ligase IRF2BP1) (EC 6.3.2.-)
Other Name:	
Human Swiss Prot No.:	Q8IU81
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IRF2BP1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year