

Antibodies Product Datasheet

Product Name:	GK5 Antibody
Catalog Number:	C49010
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GK5
Concentration:	1 mg/ml
Gene Name:	GK5
Protein Name:	Putative glycerol kinase 5 (GK 5) (Glycerokinase 5) (EC 2.7.1.30) (ATP:glycerol 3-phosphotransferase 5)
Other Name:	
Human Swiss Prot No.:	Q6ZS86
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GK5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year