

Antibodies Product Datasheet

Product Name:	FUSIP1 Antibody
Catalog Number:	C49014
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FUSIP1
Concentration:	1 mg/ml
Gene Name:	SRSF10 FUSIP1 FUSIP2 SFRS13A TASR
Protein Name:	Serine/arginine-rich splicing factor 10 (40 kDa SR-repressor protein) (SRrp40) (FUS-interacting serine-arginine-rich protein 1) (Splicing factor SRp38) (Splicing factor, arginine/serine-rich 13A) (TLS-associated protein with Ser-Arg repeats) (TASR) (TLS-a
Other Name:	
Human Swiss Prot No.:	O75494
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FUSIP1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year