

Antibodies Product Datasheet

Product Name:	GGA2 Antibody
Catalog Number:	C49015
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GGA2
Concentration:	1 mg/ml
Gene Name:	GGA2 KIAA1080
Protein Name:	ADP-ribosylation factor-binding protein GGA2 (Gamma-adaptin-related protein 2) (Golgi-localized, gamma ear-containing, ARF-binding protein 2) (VHS domain and ear domain of gamma-adaptin) (Vear)
Other Name:	
Human Swiss Prot No.:	Q9UJY4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GGA2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year