

Antibodies Product Datasheet

Product Name:	EIF2C1 Antibody
Catalog Number:	C49031
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF2C1
Concentration:	1 mg/ml
Gene Name:	EIF2C1 AGO1
Protein Name:	Protein argonaute-1 (Argonaute1) (hAgo1) (Eukaryotic translation initiation factor 2C 1) (eIF-2C 1) (eIF2C 1) (Putative RNA-binding protein Q99)
Other Name:	
Human Swiss Prot No.:	Q9UL18
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF2C1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year