

Antibodies Product Datasheet

Product Name:	DOCK3; MOCA Antibody
Catalog Number:	C49042
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DOCK3; MOCA
Concentration:	1 mg/ml
Gene Name:	DOCK3 KIAA0299 MOCA
Protein Name:	Dedicator of cytokinesis protein 3 (Modifier of cell adhesion) (Presenilin-binding protein) (PBP)
Other Name:	
Human Swiss Prot No.:	Q8IZD9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DOCK3; MOCA at Human
Storage and Stability:	Store at -20° C/1 year