

Antibodies Product Datasheet

Product Name:	CSRP2 Antibody
Catalog Number:	C49077
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CSRP2
Concentration:	1 mg/ml
Gene Name:	CSRP2 LMO5 SMLIM
Protein Name:	Cysteine and glycine-rich protein 2 (Cysteine-rich protein 2) (CRP2) (LIM domain only protein 5) (LMO-5) (Smooth muscle cell LIM protein) (SmLIM)
Other Name:	
Human Swiss Prot No.:	Q16527
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CSRP2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year