

Antibodies Product Datasheet

Product Name:	INTS9 Antibody
Catalog Number:	C49079
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human INTS9
Concentration:	1 mg/ml
Gene Name:	INTS9 RC74
Protein Name:	Integrator complex subunit 9 (Int9) (Protein related to CPSF subunits of 74 kDa) (RC-74)
Other Name:	
Human Swiss Prot No.:	Q9NV88
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of INTS9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year