

Antibodies Product Datasheet

Product Name:	H6PD Antibody
Catalog Number:	C49084
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human H6PD
Concentration:	1 mg/ml
Gene Name:	H6PD GDH
Protein Name:	GDH/6PGL endoplasmic bifunctional protein [Includes: Glucose 1-dehydrogenase (EC 1.1.1.47) (Hexose-6-phosphate dehydrogenase) 6-phosphogluconolactonase (6PGL) (EC 3.1.1.31)]
Other Name:	
Human Swiss Prot No.:	O95479
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of H6PD at Human, Mouse
Storage and Stability:	Store at -20° C/1 year