

Antibodies Product Datasheet

Product Name:	EMG1 Antibody
Catalog Number:	C49093
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EMG1
Concentration:	1 mg/ml
Gene Name:	EMG1 C2F
Protein Name:	Ribosomal RNA small subunit methyltransferase NEP1 (EC 2.1.1.-) (18S rRNA (pseudouridine-N1-)-methyltransferase NEP1) (18S rRNA Psi1248 methyltransferase) (Nucleolar protein EMG1 homolog) (Protein C2f) (Ribosome biogenesis protein NEP1)
Other Name:	
Human Swiss Prot No.:	Q92979
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EMG1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year