

Antibodies Product Datasheet

Product Name:	FUT6 Antibody
Catalog Number:	C49103
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FUT6
Concentration:	1 mg/ml
Gene Name:	FUT6 FCT3A
Protein Name:	Alpha-(1,3)-fucosyltransferase (EC 2.4.1.65) (Fucosyltransferase 6) (Fucosyltransferase VI) (Fuc-TVI) (FucT-VI) (Galactoside 3-L-fucosyltransferase)
Other Name:	
Human Swiss Prot No.:	P51993
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FUT6 at Human
Storage and Stability:	Store at -20° C/1 year