

Antibodies Product Datasheet

Product Name:	KLHL18 Antibody
Catalog Number:	C49113
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KLHL18
Concentration:	1 mg/ml
Gene Name:	KLHL18 KIAA0795 OK/SW-cl.74
Protein Name:	Kelch-like protein 18
Other Name:	
Human Swiss Prot No.:	O94889
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KLHL18 at Human
Storage and Stability:	Store at -20° C/1 year