

Antibodies Product Datasheet

Product Name:	AGA Antibody
Catalog Number:	C49128
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AGA
Concentration:	1 mg/ml
Gene Name:	AGA
Protein Name:	N(4)-(beta-N-acetylglucosaminy)-L-asparaginase (EC 3.5.1.26) (Aspartylglucosaminidase) (Glycosylasparaginase) (N4-(N-acetyl-beta-glucosaminy)-L-asparagine amidase) [Cleaved into: Glycosylasparaginase alpha chainGlycosylasparaginase beta chain]
Other Name:	
Human Swiss Prot No.:	P20933
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AGA at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year