

Antibodies Product Datasheet

Product Name:	C1QTNF3 Antibody
Catalog Number:	C49130
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1QTNF3
Concentration:	1 mg/ml
Gene Name:	C1QTNF3 CTRP3 UNQ753/PRO1484
Protein Name:	Complement C1q tumor necrosis factor-related protein 3 (Collagenous repeat-containing sequence 26 kDa protein) (CORS26) (Secretory protein CORS26)
Other Name:	
Human Swiss Prot No.:	Q9BXJ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C1QTNF3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year