

Antibodies Product Datasheet

Product Name:	C22orf41 Antibody
Catalog Number:	C49131
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C22orf41
Concentration:	1 mg/ml
Gene Name:	SYCE3 C22orf41 THEG2
Protein Name:	Synaptonemal complex central element protein 3 (Testis highly expressed gene 2 protein) (THEG-2)
Other Name:	
Human Swiss Prot No.:	A1L190
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C22orf41 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year