

Antibodies Product Datasheet

Product Name:	DNAH6 Antibody
Catalog Number:	C49132
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNAH6
Concentration:	1 mg/ml
Gene Name:	DNAH6 DNAHC6 DNHL1 HL2 KIAA1697
Protein Name:	Dynein heavy chain 6, axonemal (Axonemal beta dynein heavy chain 6) (Ciliary dynein heavy chain 6)
Other Name:	
Human Swiss Prot No.:	Q9C0G6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNAH6 at Human
Storage and Stability:	Store at -20° C/1 year