

Antibodies Product Datasheet

Product Name:	DELEC1 Antibody
Catalog Number:	C49145
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DELEC1
Concentration:	1 mg/ml
Gene Name:	DEC1 CTS9
Protein Name:	Deleted in esophageal cancer 1 (Candidate tumor suppressor CTS9)
Other Name:	
Human Swiss Prot No.:	Q9P2X7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DELEC1 at Human
Storage and Stability:	Store at -20° C/1 year