

Antibodies Product Datasheet

Product Name:	LIN7C Antibody
Catalog Number:	C49174
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LIN7C
Concentration:	1 mg/ml
Gene Name:	LIN7C MALS3 VELI3
Protein Name:	Protein lin-7 homolog C (Lin-7C) (Mammalian lin-seven protein 3) (MALS-3) (Vertebrate lin-7 homolog 3) (Veli-3)
Other Name:	
Human Swiss Prot No.:	Q9NUP9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIN7C at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year