

Antibodies Product Datasheet

Product Name:	DNASE2B Antibody
Catalog Number:	C49178
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNASE2B
Concentration:	1 mg/ml
Gene Name:	DNASE2B DLAD
Protein Name:	Deoxyribonuclease-2-beta (EC 3.1.22.1) (DNase II-like acid DNase) (DNase2-like acid DNase) (Deoxyribonuclease II beta) (DNase II beta) (Endonuclease DLAD)
Other Name:	
Human Swiss Prot No.:	Q8WZ79
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNASE2B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year