

Antibodies Product Datasheet

Product Name:	DUS3L Antibody
Catalog Number:	C49185
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DUS3L
Concentration:	1 mg/ml
Gene Name:	DUS3L
Protein Name:	tRNA-dihydrouridine(47) synthase [NAD(P)(+)]-like (EC 1.3.1.-) (tRNA-dihydrouridine synthase 3-like)
Other Name:	
Human Swiss Prot No.:	Q96G46
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DUS3L at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year