

Antibodies Product Datasheet

Product Name:	HuA/B/C/D Antibody
Catalog Number:	C49193
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HuA/B/C/D
Concentration:	1 mg/ml
Gene Name:	ELAVL4 HUD PNEM
Protein Name:	ELAV-like protein 4 (Hu-antigen D) (HuD) (Paraneoplastic encephalomyelitis antigen HuD)
Other Name:	
Human Swiss Prot No.:	P26378
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HuA/B/C/D at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year