

Antibodies Product Datasheet

Product Name:	HDDC3 Antibody
Catalog Number:	C49199
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HDDC3
Concentration:	1 mg/ml
Gene Name:	HDDC3 MESH1
Protein Name:	Guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1 (EC 3.1.7.2) (HD domain-containing protein 3) (Metazoan SpoT homolog 1) (MESH1) (Penta-phosphate guanosine-3'-pyrophosphohydrolase) ((ppGpp)ase)
Other Name:	
Human Swiss Prot No.:	Q8N4P3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HDDC3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year