

Antibodies Product Datasheet

Product Name:	HAGH Antibody
Catalog Number:	C49209
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HAGH
Concentration:	1 mg/ml
Gene Name:	HAGH GLO2 HAGH1
Protein Name:	Hydroxyacylglutathione hydrolase, mitochondrial (EC 3.1.2.6) (Glyoxalase II) (Glx II)
Other Name:	
Human Swiss Prot No.:	Q16775
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HAGH at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year