

Antibodies Product Datasheet

Product Name:	DCUN1D1 Antibody
Catalog Number:	C49213
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DCUN1D1
Concentration:	1 mg/ml
Gene Name:	DCUN1D1 DCUN1L1 RP42 SCCRO
Protein Name:	DCN1-like protein 1 (DCUN1 domain-containing protein 1) (Defective in cullin neddylation protein 1-like protein 1) (Squamous cell carcinoma-related oncogene)
Other Name:	
Human Swiss Prot No.:	Q96GG9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DCUN1D1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year