

Antibodies Product Datasheet

Product Name:	EDC3 Antibody
Catalog Number:	C49218
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EDC3
Concentration:	1 mg/ml
Gene Name:	EDC3 LSM16 YJDC YJEFN2 PP844
Protein Name:	Enhancer of mRNA-decapping protein 3 (LSM16 homolog) (YjeF N-terminal domain-containing protein 2) (YjeF_N2) (hYjeF_N2) (YjeF domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q96F86
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EDC3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year