

Antibodies Product Datasheet

Product Name:	ARD1A Antibody
Catalog Number:	C49219
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARD1A
Concentration:	1 mg/ml
Gene Name:	NAA10 ARD1 ARD1A TE2
Protein Name:	N-alpha-acetyltransferase 10 (EC 2.3.1.-) (EC 2.3.1.88) (N-terminal acetyltransferase complex ARD1 subunit homolog A) (NatA catalytic subunit)
Other Name:	
Human Swiss Prot No.:	P41227
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARD1A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year