

Antibodies Product Datasheet

Product Name:	COQ3 Antibody
Catalog Number:	C49221
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COQ3
Concentration:	1 mg/ml
Gene Name:	COQ3 UG0215E05
Protein Name:	Hexaprenyldihydroxybenzoate methyltransferase, mitochondrial (EC 2.1.1.114) (2-polyprenyl-6-hydroxyphenyl methylase) (EC 2.1.1.222) (3,4-dihydroxy-5-hexaprenylbenzoate methyltransferase) (DHMB methyltransferase) (DHMB-MT) (DHMB-MTase) (3-demethylubiquinon
Other Name:	
Human Swiss Prot No.:	Q9NZJ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COQ3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year