

Antibodies Product Datasheet

Product Name:	HADHB Antibody
Catalog Number:	C49222
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HADHB
Concentration:	1 mg/ml
Gene Name:	HADHB MSTP029
Protein Name:	Trifunctional enzyme subunit beta, mitochondrial (TP-beta) [Includes: 3-ketoacyl-CoA thiolase (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase)]
Other Name:	
Human Swiss Prot No.:	P55084
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HADHB at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year