

Antibodies Product Datasheet

Product Name:	HIBCH Antibody
Catalog Number:	C49228
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HIBCH
Concentration:	1 mg/ml
Gene Name:	HIBCH
Protein Name:	3-hydroxyisobutyryl-CoA hydrolase, mitochondrial (EC 3.1.2.4) (3-hydroxyisobutyryl-coenzyme A hydrolase) (HIB-CoA hydrolase) (HIBYL-CoA-H)
Other Name:	
Human Swiss Prot No.:	Q6NVY1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HIBCH at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year