

Antibodies Product Datasheet

Product Name:	COL4A3BP Antibody
Catalog Number:	C49229
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COL4A3BP
Concentration:	1 mg/ml
Gene Name:	COL4A3BP CERT STARD11
Protein Name:	Collagen type IV alpha-3-binding protein (Ceramide transfer protein) (hCERT) (Goodpasture antigen-binding protein) (GPBP) (START domain-containing protein 11) (StARD11) (StAR-related lipid transfer protein 11)
Other Name:	
Human Swiss Prot No.:	Q9Y5P4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COL4A3BP at Human, Mouse
Storage and Stability:	Store at -20° C/1 year