

Antibodies Product Datasheet

Product Name:	DNAJB13 Antibody
Catalog Number:	C49235
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DNAJB13
Concentration:	1 mg/ml
Gene Name:	DNAJB13 TSARG3 TSARG6
Protein Name:	DnaJ homolog subfamily B member 13 (Testis and spermatogenesis cell-related protein 6) (Testis spermatocyte apoptosis-related gene 6 protein) (Testis spermatogenesis apoptosis-related gene 3 protein) (Testis spermatogenesis apoptosis-related gene 6 protei
Other Name:	
Human Swiss Prot No.:	P59910
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DNAJB13 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year