

Antibodies Product Datasheet

Product Name:	COX4NB Antibody
Catalog Number:	C49240
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COX4NB
Concentration:	1 mg/ml
Gene Name:	EMC8 C16orf2 C16orf4 COX4AL COX4NB FAM158B NOC4
Protein Name:	ER membrane protein complex subunit 8 (Neighbor of COX4) (Protein FAM158B)
Other Name:	
Human Swiss Prot No.:	O43402
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of COX4NB at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year