

Antibodies Product Datasheet

Product Name:	KLRG1 Antibody
Catalog Number:	C49243
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KLRG1
Concentration:	1 mg/ml
Gene Name:	KLRG1 CLEC15A MAFA MAFAL
Protein Name:	Killer cell lectin-like receptor subfamily G member 1 (C-type lectin domain family 15 member A) (ITIM-containing receptor MAFA-L) (MAFA-like receptor) (Mast cell function-associated antigen)
Other Name:	
Human Swiss Prot No.:	Q96E93
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KLRG1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year