

Antibodies Product Datasheet

Product Name:	C1orf156 Antibody
Catalog Number:	C49258
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1orf156
Concentration:	1 mg/ml
Gene Name:	METTL18 ASTP2 C1orf156
Protein Name:	Histidine protein methyltransferase 1 homolog (EC 2.1.1.-) (Arsenic-transactivated protein 2) (AsTP2) (Methyltransferase-like protein 18)
Other Name:	
Human Swiss Prot No.:	O95568
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C1orf156 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year