

Antibodies Product Datasheet

Product Name:	Band 3/AE1 Antibody
Catalog Number:	C49282
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Band 3/AE1
Concentration:	1 mg/ml
Gene Name:	SLC4A1 AE1 DI EPB3
Protein Name:	Band 3 anion transport protein (Anion exchange protein 1) (AE 1) (Anion exchanger 1) (Solute carrier family 4 member 1) (CD antigen CD233)
Other Name:	
Human Swiss Prot No.:	P02730
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Band 3/AE1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year