

Antibodies Product Datasheet

Product Name:	FBXO11 Antibody
Catalog Number:	C49294
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein. AA range:1500-1500
Concentration:	1 mg/ml
Gene Name:	FBXO11 FBX11 VIT1 UG063H01
Protein Name:	F-box only protein 11 (Vitiligo-associated protein 1) (VIT-1)
Other Name:	
Human Swiss Prot No.:	Q86XK2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FBXO11 at Human,Rat
Storage and Stability:	Store at -20° C/1 year