

Antibodies Product Datasheet

Product Name:	Alkyl-DHAP synthase/AGPS Antibody
Catalog Number:	C49304
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Alkyl-DHAP synthase/AGPS
Concentration:	1 mg/ml
Gene Name:	AGPS AAG5
Protein Name:	Alkyldihydroxyacetonephosphate synthase, peroxisomal (Alkyl-DHAP synthase) (EC 2.5.1.26) (Aging-associated gene 5 protein) (Alkylglycerone-phosphate synthase)
Other Name:	
Human Swiss Prot No.:	O00116
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Alkyl-DHAP synthase/AGPS at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year