

Antibodies Product Datasheet

Product Name:	C1orf77 Antibody
Catalog Number:	C49306
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1orf77
Concentration:	1 mg/ml
Gene Name:	CHTOP C1orf77 FOP HT031 PP7704
Protein Name:	Chromatin target of PRMT1 protein (Friend of PRMT1 protein) (Small arginine- and glycine-rich protein) (SRAG)
Other Name:	
Human Swiss Prot No.:	Q9Y3Y2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C1orf77 at Human
Storage and Stability:	Store at -20° C/1 year