

Antibodies Product Datasheet

Product Name:	ALKBH1 Antibody
Catalog Number:	C49312
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ALKBH1
Concentration:	1 mg/ml
Gene Name:	ALKBH1 ABH ABH1 ALKBH
Protein Name:	Alkylated DNA repair protein alkB homolog 1 (EC 1.14.11.33) (Alpha-ketoglutarate-dependent dioxygenase ABH1) (DNA lyase ABH1) (EC 4.2.99.18) (DNA oxidative demethylase ALKBH1)
Other Name:	
Human Swiss Prot No.:	Q13686
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ALKBH1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year