

Antibodies Product Datasheet

Product Name:	AP3D1 Antibody
Catalog Number:	C49314
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AP3D1
Concentration:	1 mg/ml
Gene Name:	AP3D1 PRO0039
Protein Name:	AP-3 complex subunit delta-1 (AP-3 complex subunit delta) (Adapter-related protein complex 3 subunit delta-1) (Delta-adaptin)
Other Name:	
Human Swiss Prot No.:	O14617
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AP3D1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year