

Antibodies Product Datasheet

Product Name:	EPB41L1 Antibody
Catalog Number:	C49321
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EPB41L1
Concentration:	1 mg/ml
Gene Name:	EPB41L1 KIAA0338
Protein Name:	Band 4.1-like protein 1 (Neuronal protein 4.1) (4.1N)
Other Name:	
Human Swiss Prot No.:	Q9H4G0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EPB41L1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year