

Antibodies Product Datasheet

Product Name:	KCNIP2 Antibody
Catalog Number:	C49325
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KCNIP2
Concentration:	1 mg/ml
Gene Name:	KCNIP2 KCHIP2
Protein Name:	Kv channel-interacting protein 2 (KChIP2) (A-type potassium channel modulatory protein 2) (Cardiac voltage-gated potassium channel modulatory subunit) (Potassium channel-interacting protein 2)
Other Name:	
Human Swiss Prot No.:	Q9NS61
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KCNIP2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year