

Antibodies Product Datasheet

Product Name:	KCNQ3 Antibody
Catalog Number:	C49326
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KCNQ3
Concentration:	1 mg/ml
Gene Name:	KCNQ3
Protein Name:	Potassium voltage-gated channel subfamily KQT member 3 (KQT-like 3) (Potassium channel subunit alpha KvLQT3) (Voltage-gated potassium channel subunit Kv7.3)
Other Name:	
Human Swiss Prot No.:	O43525
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KCNQ3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year