

Antibodies Product Datasheet

Product Name:	KMT2C Antibody
Catalog Number:	C49327
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KMT2C
Concentration:	1 mg/ml
Gene Name:	MLL3 HALR KIAA1506 KMT2C
Protein Name:	Histone-lysine N-methyltransferase MLL3 (EC 2.1.1.43) (Homologous to ALR protein) (Lysine N-methyltransferase 2C) (KMT2C) (Myeloid/lymphoid or mixed-lineage leukemia protein 3)
Other Name:	
Human Swiss Prot No.:	Q8NEZ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KMT2C at Human
Storage and Stability:	Store at -20° C/1 year