

Antibodies Product Datasheet

Product Name:	C6orf182 Antibody
Catalog Number:	C49357
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C6orf182
Concentration:	1 mg/ml
Gene Name:	CEP57L1 C6orf182 CEP57R
Protein Name:	Centrosomal protein CEP57L1 (Centrosomal protein 57kDa-like protein 1) (Centrosomal protein of 57 kDa-related protein) (Cep57R) (Cep57-related protein)
Other Name:	
Human Swiss Prot No.:	Q8IYX8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C6orf182 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year