

Antibodies Product Datasheet

Product Name:	B3GALNT2 Antibody
Catalog Number:	C49365
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human B3GALNT2
Concentration:	1 mg/ml
Gene Name:	B3GALNT2
Protein Name:	UDP-GalNAc:beta-1,3-N-acetylgalactosaminyltransferase 2 (Beta-1,3-GalNAc-T2) (EC 2.4.1.-) (Beta-1,3-N-acetylgalactosaminyltransferase II)
Other Name:	
Human Swiss Prot No.:	Q8NCR0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of B3GALNT2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year