

Antibodies Product Datasheet

Product Name:	IL-17RE Antibody
Catalog Number:	C49373
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IL-17RE
Concentration:	1 mg/ml
Gene Name:	IL17RE UNQ3056/PRO9877
Protein Name:	Interleukin-17 receptor E (IL-17 receptor E) (IL-17RE)
Other Name:	
Human Swiss Prot No.:	Q8NFR9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IL-17RE at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year