

Antibodies Product Datasheet

Product Name:	CCDC68 Antibody
Catalog Number:	C49406
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCDC68
Concentration:	1 mg/ml
Gene Name:	CCDC68
Protein Name:	Coiled-coil domain-containing protein 68 (Cutaneous T-cell lymphoma-associated antigen se57-1) (CTCL-associated antigen se57-1)
Other Name:	
Human Swiss Prot No.:	Q9H2F9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCDC68 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year